"I was bitten by a wasp on the job site. It was really painful but I used the bite away and the pain was gone in seconds and only an hour later when I checked the

swelling had gone as well. Now I always have one in my work truck."

Ashley, Electrician.

"My 9 yr old granddaughter came up to me and said "can I have that thing that stops the itch" I was so pleased that we have found something that she will use and it stop the itching."

Keryn, Editor outbacktravelaustralia.com.au



"My child comes up in red welts when bitten by Mosquitoes so we have to be so careful. Now we use the bite away and no more reaction,"

Cindy, Mum.

"I recently used the bite away after being bitten by sand flies. Usually I react really badly with redness and itching for days. One application to each bite and the itching stopped. It was quick,

painless and easy to use. I highly recommend bite away as a useful addition to the home and travel medical kit."

> Libby, Occupational therapist.



Relief in seconds

Medical device for the external treatment of insect bites & stings.

Bite Away works on natural concentrated heat being applied directly to the bite site.

An integrated microprocessor allow the temperature to reach approx. 50°C and remain constant for 3 or 6 seconds.

Heat being applied to the bite denatures the toxins/proteins and therefore stops the histomine reaction.

NO SWELLING, NO PAIN, NO ITCHING

As a result allergic reactions to bites are prevented and bites heal very quickly.

Works on 2 x AA Batteries which equals 300 applications. Then just change the batteries and keep healing.

Exclusive importers Bite Away Australia, P.O. Box 1024, Byron Bay, NSW, 2481 1300 559 519 www.biteawayaustralia.com.au

German quality manufactured by

RIEMSER 4





1300 559 519

biteawayaustralia.com.au









treats...

- midge/sand fly bites
- mosquitoes bites
- wasp stings
- bee stings
- flea bites
- lice
- bed bugs
- stinging nettle
- blue bottle jelly fish stings
- horse fly/march fly bites
- mites/chigger bites
- hairy grub

helps with...

- tick bites
- green ants stings
- any insect bite with a histamine reaction





MIDGE BITES

IN ONLY 6

SECONDS









250,000 sold worldwide

- relieves swelling, pain and itching
- prevents allergic reaction
- portable, easy to use
- patented technology
- clinically tested
- heals bites and stings in seconds
- chemical free







safe for...

- all your family and pets
- people with heart conditions, epilepsy or pacemakers
- people with allergies
- pregnant women

Instructions...

- Select button of your choice.
- Place round contact plate on bite.
- Press button.
- Wait for LED light treatment has started.
- When LED light goes out treatment has ended.
- Remove after each treatment.
- Reapply if necessary.

Do not use on mucous membranes, open wounds or eyes.

Stops itching safely and fast at the press of a button.

A chemical free alternative for the instant relief of insect bites & stings.









